



SKU:

SUNRISE MEDICAL BREEZY RUBIX WHEELCHAIR

\$1,999.99

The Sunrise Medical Breezy Rubix Wheelchair is a folding, heavy-duty wheelchair that supports individuals up to 170 kg. With seat widths up to 24" and multiple adjustable features, it offers tailored comfort, postural support and reliable daily mobility.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Sizes
S10000	\$1,999.99	The Sunrise Medical Breezy Rubix Wheelchair is a folding, heavy-duty wheelchair that supports individuals up to 170 kg. With seat widths up to 24" and multiple adjustable features, it offers tailored comfort, postural support and reliable daily mobility.	20"
S10001	\$1,999.99	The Sunrise Medical Breezy Rubix Wheelchair is a folding, heavy-duty wheelchair that supports individuals up to 170 kg. With seat widths up to 24" and multiple adjustable features, it offers tailored comfort, postural support and reliable daily mobility.	22"
S10002	\$1,999.99	The Sunrise Medical Breezy Rubix Wheelchair is a folding, heavy-duty wheelchair that supports individuals up to 170 kg. With seat widths up to 24" and multiple adjustable features, it offers tailored comfort, postural support and reliable daily mobility.	24"

PRODUCT DESCRIPTION

The Breezy Rubix Wheelchair is designed for individuals needing wider seating and a higher weight capacity. It features a reinforced aluminium frame, adjustable seat depth, configurable seat height, and a tension-adjustable backrest for tailored support.

Quick-release rear wheels and swing-away, angle-adjustable legrests support easier transfers and compact storage.

Additional features include flip-back armrests, anti-tip wheels and rear wheel park brakes for everyday usability and safety. It also has crash test certification.

Who it's for

- For clients seeking a reliable, adjustable heavy-duty wheelchair with wider seat width and self-propulsion.
- For therapists prescribing a bariatric wheelchair with certified safety, daily durability and adjustability.

Key features

- **Adjustable seat depth (46 cm - 51 cm)** – supports postural alignment and comfort
 - **Heavy-duty aluminium frame** – accommodates up to 170 kg
 - **Available seat widths** – 20", 22" and 24" accommodates larger body shapes
 - **Flip-back armrests** – allow safer lateral transfers and access
 - **Tension-adjustable backrest** – offers personalised support
 - **Swing-away, angle-adjustable legrests** – enhance positioning and transfers
 - **Quick-release rear wheels** – for ease of transport and handling
 - **Rear park brake** – accessible and easy to operate
 - **Crash tested** – certified for transport safety

ADDITIONAL INFORMATION

Sizes	20", 22", 24"
Safe Working Load (SWL)	170 kg
Net Weight	16.8 kg
Seat Width	20" / 22" / 24"
Seat Depth	46 cm - 51 cm
Seat Height (Floor to Seat)	37.5 cm - 51 cm
Backrest Height (From Seat)	43.5 cm (adjustable 41 cm - 46cm)
Armrest Height (From Seat)	20 cm
Front Castor	8"
Rear Wheel	24"
Brand	Sunrise Medical

[Breezy Rubix User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

No FAQ Found