



SCALEO PREMIUM PATIENT LIFTER

\$8,995.99 - \$10,295.99

5 years on frame, 2 yrs on electronics, 1 year on battery

Features:

- 4-point overhead suspension or an electrical tilting suspension
- SLS hook - Safe Locking System
- Removable high capacity battery
 - Electrical opening of the base
 - Lift patient from floor
- Mechanical and electrical emergency lowering
- Lightweight aluminium structure



SKU: S13721

Categories: [Floor Hoist and Lifters](#), [Patient Transfer Equipment](#)

VARIATIONS

Image

SKU

Price

Description

Weight Scale



S13721

~~\$10,995.99~~ \$10,295.99

With Weight Scale



S16989

~~\$9,995.99~~ \$8,995.99

Without Weight Scale

GALLERY IMAGES



PRODUCT DESCRIPTION

New model available in autumn Poweo® 200 is dedicated to all kind of transfer needs in a hospital environment, from a daily transfer up to the most specialized transfer in Intensive Care Units.

It can be equipped with a 4-point overhead suspension or an electrical tilting suspension allowing transfer without any effort for the caregiver. Poweo 200 can be used with various accessories such as a scale, a stretcher or a large range of slings adapted to all kind of morphologies and pathologies.

ADDITIONAL INFORMATION

Weighing capacity	200 kg
NiMh battery capacity	> 90 cycles
Actuator	24 V - 10 A
NiMh Battery	24 V - 3 Ah
Front wheels	ø 75 mm
Rear wheels with brakes	ø 100 mm
Width legs closed	68 cm
Width legs open	118 cm
Length in use	132 cm
Maximum height folded	145 cm
Maximum height in use	214 cm
Height under suspension / with Scalis	187 cm / 175 cm (4 points suspension), 160 cm / 147 cm (Electric Suspensions)
Legs height	11 cm

Ground clearance	5 cm
Total weight	49 kg (4 points suspension), 64 kg (electric suspension)
Material	Aluminium and steel
Brand	Scaleo