



SKU:

DRIVE MEDICAL SCOUT MOBILITY SCOOTER

\$2,199.99

The Drive Medical Scout Mobility Scooter is a lightweight, transportable scooter designed for comfort and convenience. With a 20 km range and tool-free dismantling, it supports everyday mobility needs.



Scan QR code to
view product online

VARIATIONS

Image

SKU

Price

Description

Product Colour



S11712 \$2,199.99

Red



S11808 \$2,199.99

Blue

PRODUCT DESCRIPTION

The Drive Medical Scout Mobility Scooter is a compact four-wheel scooter designed for portability without compromising on comfort or functionality. Ideal for everyday use, it easily disassembles into five lightweight components – seat, front and rear sections, battery pack and basket – making it a practical choice for transport and storage.

A 360° swivelling seat, width-adjustable armrests and an adjustable tiller offer a customisable setup to enhance comfort during use. The lightweight, padded seat adjusts in height without tools, supporting ergonomic positioning.

The Scout's removable battery pack includes a carry handle and lock, allowing both on-board and off-board charging. With a travel range of up to 20 km and top speed of 6.5 km/h, it provides dependable performance for daily errands and short trips.

Safety features include rear anti-tip wheels, an automatic electromagnetic braking system, and non-marking puncture-proof tyres.

Who it's for

- For clients seeking a reliable, lightweight scooter that's easy to transport and store.
- For therapists recommending a value-focused mobility solution with practical features and safety controls.

Key features

- **Splits into 5 lightweight components** – simplifies transport and storage
- **Up to 20 km travel range** – ideal for daily short trips
- **360° swivelling padded seat** – supports safe and easy transfers
- **Height-adjustable seat and adjustable tiller** – allows personal comfort customisation
- **Delta handlebars and width-adjustable armrests** – enhances control and ergonomic use
- **Removable battery with carry handle** – offers charging flexibility
- **Electromagnetic braking system** – increases rider safety
- **Anti-tip wheels and non-marking tyres** – for stable and safe travel across surfaces
- **Front storage basket** – provides a convenient way to carry essentials

ADDITIONAL INFORMATION

Product Safe Working Load (SWL) - in kg	130 kg
Product Battery Size / Type - in AH / (gel/lead acid/lithium)	2 x 20 AH Gel
	Product Range - in km 20 km
	Product Speed (Max) - in km/hr 6 km/hr
Product Range - in km	20 km
Product Speed (Max) - in km/hr	6 km/hr
	Mechanical 12 months
Product Warranty - in months	Electrical 12 months
	Batteries 12 months"
Product Size - in carton(s) (WxLxH) in cm	50 x 133 x 51
Product Size - in carton(s) (gross weight) in kg	63 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	48 x 108 x 91
Product Size - overall assembled (gross weight) in kg	48 kg
Product Size - seat base width in cm	45 cm
	Mechanical 12 months
Product Warranty	Electrical 12 months
	Batteries 12 months

Product Size - seat base depth in cm	44 cm
Product Size - seat back width in cm	44 cm
Product Size - seat back height in cm	35.5 cm
Product Size - seat to floor in cm	43 cm – 48 cm
Product Size - front/ rear tyre type/size	8" / Solid
Product Size - front/rear wheel size in inches	5"
Product Drive - turning circle in cm	140 cm
Product Drive - max safe ground clearance in cm	6 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	No
Product Charger Type	1.8 AH
Product Motor Size	2 Pole / 270 Watt
Product Display	Analog