

SKU:

DRIVE MEDICAL SCOUT PORTABLE MOBILITY SCOOTER

\$1,899.99

The Drive Medical Scout portable mobility scooter is big on features including the next generation Drive Medical splitting mechanism for easy storage and transportation. With a maximum distance range of 16km on a full battery charge, it has been developed for superb comfort, style, performance, and reliability.



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Product Colour
	S11712	\$1,899.99		Red
	S11808	\$1,899.99		Blue

PRODUCT DESCRIPTION

Drive Medical Scout Portable Mobility Scooter is an amazing and portable scooter able to move you anywhere and folds to be tucked away . Offering superior value with many easy-to-use functions, the Drive Medical Scout mobility scooter is lightweight and stylish. It is compact in size and splits easily for storage and transportation.

Key features include:

- Maximum range of 16km on full battery charge
- Comfortable lightweight seat which swivels 360°
- Removable battery pack with easy to carry handle and locking mechanism
- Top speed of 6.5kph (4mph)
- On board and in-house charging as standard
- Automatic electromagnetic brake system
- Tool free height adjustable seat which is easy to remove
- Delta bars as standard
- Easy-to-change body panels
- Adjustable tiller which can fold down
- Stylish silver wheel hubs
- Black puncture proof non-marking tyres
- Width adjustable armrests
- Front removable storage basket as standard
- Rear anti tip wheels as standard.

ADDITIONAL INFORMATION

Product Safe Working Load (SWL) - in kg	130 kg
Product Battery Size / Type - in AH / (gel/lead acid/lithium)	2 x 20 AH Gel
	Product Range - in km 20 km
	Product Speed (Max) - in km/hr 6 km/hr
Product Range - in km	20 km
Product Speed (Max) - in km/hr	6 km/hr
	Mechanical 12 months
Product Warranty - in months	Electrical 12 months
	Batteries 12 months"
Product Size - in carton(s) (WxLxH) in cm	50 x 133 x 51
Product Size - in carton(s) (gross weight) in kg	63 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	48 x 108 x 91
Product Size - overall assembled (gross weight) in kg	48 kg
Product Size - seat base width in cm	45 cm
	Mechanical 12 months
Product Warranty	Electrical 12 months
	Batteries 12 months

Product Size - seat base depth in cm	44 cm
Product Size - seat back width in cm	44 cm
Product Size - seat back height in cm	35.5 cm
Product Size - seat to floor in cm	43 cm – 48 cm
Product Size - front/ rear tyre type/size	8" / Solid
Product Size - front/rear wheel size in inches	5"
Product Drive - turning circle in cm	140 cm
Product Drive - max safe ground clearance in cm	6 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	No
Product Charger Type	1.8 AH
Product Motor Size	2 Pole / 270 Watt
Product Display	Analog