


SKU:

## **DRIVE MEDICAL SCOUT PORTABLE MOBILITY SCOOTER**

**\$1,899.99**



The Drive Medical Scout portable mobility scooter is big on features including the next generation Drive Medical splitting mechanism for easy storage and transportation. With a maximum distance range of 16km on a full battery charge, it has been developed for superb comfort, style, performance, and reliability.



Scan QR code to  
view product online



VARIATIONS

Image	SKU	Price	Description	Product Colour
	S11712	\$1,899.99		Red
	S11808	\$1,899.99		Blue

## **PRODUCT DESCRIPTION**

Drive Medical Scout Portable Mobility Scooter is an amazing and portable scooter able to move you anywhere and folds to be tucked away . Offering superior value with many easy-to-use functions, the Drive Medical Scout mobility scooter is lightweight and stylish. It is compact in size and splits easily for storage and transportation.

Key features include:

- Maximum range of 16km on full battery charge
- Comfortable lightweight seat which swivels 360°
- Removable battery pack with easy to carry handle and locking mechanism
- Top speed of 6.5kph (4mph)
- On board and in-house charging as standard
- Automatic electromagnetic brake system
- Tool free height adjustable seat which is easy to remove
- Delta bars as standard
- Easy-to-change body panels
- Adjustable tiller which can fold down
- Stylish silver wheel hubs
- Black puncture proof non-marking tyres
- Width adjustable armrests
- Front removable storage basket as standard
- Rear anti tip wheels as standard.

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL) - in kg</b>	130 kg
<b>Product Battery Size / Type - in AH / (gel/lead acid/lithium)</b>	2 x 20 AH Gel
	Product Range - in km 20 km
	Product Speed (Max) - in km/hr 6 km/hr
<b>Product Range - in km</b>	20 km
<b>Product Speed (Max) - in km/hr</b>	6 km/hr
<b>Product Warranty - in months</b>	Mechanical 12 months Electrical 12 months Batteries 12 months"
<b>Product Size - in carton(s) (WxLxH) in cm</b>	50 x 133 x 51
<b>Product Size - in carton(s) (gross weight) in kg</b>	63 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	48 x 108 x 91
<b>Product Size - overall assembled (gross weight) in kg</b>	48 kg
<b>Product Size - seat base width in cm</b>	45 cm
<b>Product Warranty</b>	Mechanical 12 months Electrical 12 months Batteries 12 months

<b>Product Size - seat base depth in cm</b>	44 cm
<b>Product Size - seat back width in cm</b>	44 cm
<b>Product Size - seat back height in cm</b>	35.5 cm
<b>Product Size - seat to floor in cm</b>	43 cm – 48 cm
<b>Product Size - front/ rear tyre type/size</b>	8" / Solid
<b>Product Size - front/rear wheel size in inches</b>	5"
<b>Product Drive - turning circle in cm</b>	140 cm
<b>Product Drive - max safe ground clearance in cm</b>	6 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	No
<b>Product Charger Type</b>	1.8 AH
<b>Product Motor Size</b>	2 Pole / 270 Watt
<b>Product Display</b>	Analog