



**SKU:**

## **DRIVE MEDICAL SCOUT MOBILITY SCOOTER**

**\$1,999.99**

The Drive Medical Scout Mobility Scooter is a lightweight, transportable scooter designed for comfort and convenience. With a 20 km range and tool-free dismantling, it supports everyday mobility needs.



Scan QR code to  
view product online

## VARIATIONS

**Image**

**SKU**

**Price**

**Description**

**Product Colour**



S11712 \$1,999.99

Red



S11808 \$1,999.99

Blue

## PRODUCT DESCRIPTION

The Drive Medical Scout Mobility Scooter is a compact four-wheel scooter designed for portability without compromising on comfort or functionality. Ideal for everyday use, it easily disassembles into five lightweight components – seat, front and rear sections, battery pack and basket – making it a practical choice for transport and storage.

A 360° swivelling seat, width-adjustable armrests and an adjustable tiller offer a customisable setup to enhance comfort during use. The lightweight, padded seat adjusts in height without tools, supporting ergonomic positioning.

The Scout's removable battery pack includes a carry handle and lock, allowing both on-board and off-board charging. With a travel range of up to 20 km and top speed of 6.5 km/h, it provides dependable performance for daily errands and short trips.

Safety features include rear anti-tip wheels, an automatic electromagnetic braking system, and non-marking puncture-proof tyres.

### Who it's for

- For clients seeking a reliable, lightweight scooter that's easy to transport and store.
- For therapists recommending a value-focused mobility solution with practical features and safety controls.

### Key features

- **Splits into 5 lightweight components** – simplifies transport and storage
- **Up to 20 km travel range** – ideal for daily short trips
- **360° swivelling padded seat** – supports safe and easy transfers
- **Height-adjustable seat and adjustable tiller** – allows personal comfort customisation
- **Delta handlebars and width-adjustable armrests** – enhances control and ergonomic use
- **Removable battery with carry handle** – offers charging flexibility
- **Electromagnetic braking system** – increases rider safety
- **Anti-tip wheels and non-marking tyres** – for stable and safe travel across surfaces
- **Front storage basket** – provides a convenient way to carry essentials

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL) - in kg</b>	130 kg
<b>Product Battery Size / Type - in AH / (gel/lead acid/lithium)</b>	2 x 20 AH Gel
	Product Range - in km 20 km
	Product Speed (Max) - in km/hr 6 km/hr
<b>Product Range - in km</b>	20 km
<b>Product Speed (Max) - in km/hr</b>	6 km/hr
<b>Product Warranty - in months</b>	Mechanical 12 months Electrical 12 months Batteries 12 months"
<b>Product Size - in carton(s) (WxLxH) in cm</b>	50 x 133 x 51
<b>Product Size - in carton(s) (gross weight) in kg</b>	63 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	48 x 108 x 91
<b>Product Size - overall assembled (gross weight) in kg</b>	48 kg
<b>Product Size - seat base width in cm</b>	45 cm
<b>Product Warranty</b>	Mechanical 12 months Electrical 12 months Batteries 12 months

<b>Product Size - seat base depth in cm</b>	44 cm
<b>Product Size - seat back width in cm</b>	44 cm
<b>Product Size - seat back height in cm</b>	35.5 cm
<b>Product Size - seat to floor in cm</b>	43 cm – 48 cm
<b>Product Size - front/ rear tyre type/size</b>	8" / Solid
<b>Product Size - front/rear wheel size in inches</b>	5"
<b>Product Drive - turning circle in cm</b>	140 cm
<b>Product Drive - max safe ground clearance in cm</b>	6 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	No
<b>Product Charger Type</b>	1.8 AH
<b>Product Motor Size</b>	2 Pole / 270 Watt
<b>Product Display</b>	Analog