

SKU: S13390



DRIVE MEDICAL SCOUT SPORT QUATTRO MOBILITY SCOOTER

From: \$2,199.99

The Drive Medical Scout Sport Quattro is a fantastic deluxe portable mobility scooter for people who have active lifestyles and enjoy travel. It has superior features such as all-round suspension and the next generation Drive Medical splitting mechanism for easy storage and transportation.

With a maximum distance range of 20km on a full battery charge, it is suited to people who are looking for a comfortable, high-performance mobility scooter that can be used anywhere.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Drive Medical - Scout Sport Quattro Mobility Scooter

The Scout Sport Quattro is a stylish and sporty portable mobility scooter. Its next generation splitting mechanism allows the seat to be removed and the tiller to fold down so that you can fit it in the boot of the car, or take it for trips on trains, cruises, or airplanes.

It is easy to operate with many features to ensure a comfortable and safe ride. In particular, the stylish, black alloy wheels have puncture proof tyres for increased durability.

Key features include:

- Maximum range of 20km on full battery charge
- All-round suspension – The stylish silver wheel hubs with black puncture proof non-marking tyres ensures a comfortable ride for a travelling range of up to 20 kilometres on a full battery.
- Top speed of 6.5km/hr
- Next generation splitting mechanism includes removable seat and fold down tiller for easy storage and transportation
- Superior safety features - LED front light for maximum visibility, automatic electromagnetic brake system, anti-roll back system safety device which prevents the mobility scooter from rolling backwards on hills, rear anti-tip wheels
- Comfortable lightweight seat – deluxe cushioned seating with stitching which swivels 360° for easy transfers on and off the scooter. The seat is also height adjustable to suit your comfort needs.
- Freewheel facility allows movement of the mobility scooter without switching the motor on
- On board and in house charging
- Delta handlebars for easier use
- Adjustable tiller
- Flip-up contoured armrests which are width adjustable
- Removable battery pack – fitted with airline friendly batteries, it has an easy-to-carry handle and easy-to-use locking mechanism
- Front removable storage basket

ADDITIONAL INFORMATION

Brand	Drive Medical
Range	Up to 20km
Maximum speed	6.5km/hr
Safe working load (SWL)	130kg (286.6 lbs)
Overall assembled dimensions	48cm (W) x 106cm (L) x 91cm (H)
Overall assembled gross weight	47.7kg
Seat width	45cm
Seat depth	44cm
Seat back - height	35.5cm
Seat to floor - height	54cm – 59cm
Front/rear tyre type/size	8" / Solid
Front/rear wheel size	5" inches
Drive - turning circle in cm	140cm
Maximum ground clearance	6cm
Maximum incline slope	8 degrees
Drive - suspension	Yes (front and rear)
Battery type	2 x 20 AH gel batteries
Charger type	1.8 AH
Motor size	2 Pole / 270 watt

Product lighting	LED
Product display	Analog
Colour	Red
Warranty	12 months for mechanical and electrical components 12 months for batteries