

SKU: S13390

DRIVE MEDICAL SCOUT SPORT QUATTRO MOBILITY SCOOTER

\$2,499.99

Exclusive to ILS, the Drive Medical Scout Sport Quattro Mobility Scooter offers advanced comfort and portability in a compact design. With all-round suspension and tool-free dismantling, it supports active lifestyles and smoother rides.



Scan QR code to
view product online

GALLERY IMAGES





PRODUCT DESCRIPTION

The Drive Medical Scout Sport Quattro Mobility Scooter is a deluxe portable scooter built for comfort, safety and flexibility. Designed for those with active routines, it features all-round suspension to absorb vibration and enhance the ride experience across uneven surfaces.

Its next-generation splitting mechanism enables quick disassembly into five manageable components – seat, front and rear sections, battery pack and basket – making it easy to transport in a car boot or store when not in use.

Comfort features include a deluxe padded seat with 360° swivel, adjustable tiller and width-adjustable armrests with contoured padding. The seat height can also be changed without tools.

Additional features include a bright LED headlight for visibility, electromagnetic brakes with an anti-rollback function, and a removable battery pack with charging flexibility.

Compact yet capable, the Scout Sport Quattro offers a top speed of 6.5 km/h and a range of up to 20 km on a full charge.

Who it's for

- For clients seeking a compact scooter with full suspension and easy transport features.
- For therapists wanting a portable option with enhanced ride comfort and safety components.

Key features

- **Splits into 5 portable components** – for easy transport and storage
- **All-round suspension** – delivers a smoother, more comfortable ride
 - **360° swivelling deluxe seat** – supports easier, safer transfers
 - **Adjustable tiller and seat height** – allows individualised setup
- **Width-adjustable padded armrests** – improve positioning and support
- **LED headlight and electromagnetic brake system** – promotes safety and visibility
- **Removable battery pack with locking handle** – provides flexible charging options

- **Delta handlebars and anti-tip rear wheels** – for controlled, stable handling

ADDITIONAL INFORMATION

Safe Working Load (SWL)	130 kg
Battery Type	Sealed Lead Acid
Battery Size	2 x 22 AH
Range	20 km
Max Speed	6.5 km/h
Total Weight	Seat - 7 kg Front Section - 16 kg Rear Section - 10 kg Battery Pack - 15 kg Total Weight - 47 kg
Suspension	Front & Rear
Display	Analog
Max Safe Incline Slope	8°
Turning Circle	140 cm
Ground Clearance	64 cm
Brand	Drive Medical

[Drive Medical Scout Sport Quattro User Manual](#)

[Drive Medical Scout Sport Quattro MSDS Certificate](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

WHAT ARE THE SIZE OF THE WHEELS AND ARE THEY SOLID AND WHAT GROUND CLEARANCE OF THE SCOUT STYLE SCOOTER?

The Scout Style Mobility Scooter has solid tires with the sizing of 20 x 50cm (8 x 2") . The Ground Clearance is 6cm (2.5")

WE PLAN ON TRAVELLING A LOT. WHAT TYPE OF BATTERIES DOES THE SCOUT STYLE PLUS SCOOTER HAVE?

The Scout Style Mobility Scooter is equipped with **2 x 20AH airline friendly** batteries.

I NEED TO KNOW IF THE SCOOTER WILL FIT INTO MY CAR. CAN I HAVE THE DIMENSIONS OF THE SCOUT STYLE MOBILITY SCOOTER?

The Scout Style Mobility Scooter dimensions are as follows,

Length	106cm
Height	91cm
Width	48cm

COULD YOU LET ME KNOW THE MAX USER WEIGHT FOR THE SCOUT STYLE MOBILITY SCOOTER?

The maximum user weight for the Scout Style Mobility Scooter is 130kg

WILL I BE ABLE TO EASILY LIFT THE SCOUT STYLE MOBILITY SCOOTER FROM MY VEHICLE? HOW MUCH DOES IT ACTUALLY WEIGH?

The Scout Style Mobility Scooter weighs 45kg. You may need assistance to lift the scooter.