


SKU: S0234

SEAT WALKER - HEAVY DUTY HIGH MACK

\$399.99

Seat Walker constructed from heavy-duty aluminium for durability and features a removable padded backrest with easy release pins. Handles are ergonomically designed and the handbrakes can be locked for stability. The padded seat allows for comfort and the shopping basket is a handy aid and can be removed when folding for transportation and storage. The 8" x 2" castors enable easy rolling over varying terrain.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Heavy Duty High Mack was designed for taller people in mind, it's handle height range can go from 870 mm – 980 mm and can take an impressive 150kg in weight. (please look at the specifications tab for more information)

ADDITIONAL INFORMATION

Brand	Ausmedic
Max SWL	150 kg (330 lbs)
Overall width	70 cm (27.5 inches)
Overall depth	64 cm (25.1 inches)
Floor to seat height	58 cm (22.8 inches)
Seat width	50 cm (19.6 inches)
Handle height range	87 cm – 98 cm 34.2 – 38.5 inches
Weight of walker	9.4 kg (20.7 lbs)
Colour	Red
Warranty	12 Months Warranty. Please note: Brake maintenance and adjustment is not covered under warranty.