

SKU:

HERO MEDICAL 6" STANDARD SEAT WALKER

\$169.99

The Hero Medical 6" Standard Seat Walker is a compact, four-wheeled walking aid that offers safe, stable mobility support indoors and outdoors. With a padded seat, height-adjustable handles and built-in storage, it supports everyday independence.



Scan QR code to
view product online



VARIATIONS






Image	SKU	Price	Description	Colour
	S11167	\$169.99		Burgundy
	S11991	\$169.99		Black
	S11990	\$169.99		Champagne
	S11168	\$169.99		Blue

Image	SKU	Price	Description	Colour
	S11169	\$169.99		Purple

PRODUCT DESCRIPTION

The Hero Medical 6" Standard Seat Walker is a practical mobility aid for individuals needing assistance with walking and resting. Its lightweight aluminium frame folds easily for transport and is fitted with 6" wheels to navigate a variety of surfaces, whether indoors or outdoors.

Height-adjustable handles help ensure ergonomic positioning, while light-touch loop brakes enhance control and stopping ease.

A cushioned seat and backrest support comfort during breaks, reducing fatigue and promoting safety.

An under-seat storage bag provides convenient space for carrying personal belongings, making this walker a functional solution for daily use.

Who it's for

- For clients who need a cost-effective, standard walker to support mobility and reduce fatigue.
- For therapists prescribing a durable, compact walker with adjustable handbrakes to assist safe gait patterns.

Key features

- **Lightweight aluminium frame** – easy to manoeuvre and store
- **Four-wheeled design** – offers stability and smooth movement
- **Padded seat and backrest** – support rest and comfort during use
- **Adjustable handle height** – promotes correct posture and safe walking
- **Light-touch loop brakes** – improve safety and control
- **Folding frame with storage bag** – allows easy transport and keeps items accessible

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight capacity / Max SWL	140kg (308 lbs)
Overall height (floor to handles)	77 cm – 92 cm
Overall length	69 cm (27 inches)
Overall width	62 cm (24.4 inches)
Seat height	55 cm (21.6 inches)
Seating width	44 cm (17.3 inches)
Seat depth	33 cm (12.9 inches)
Folded length	82cm (32.2 inches)
Folded width	23 cm (9 inches)
Outer width (including back castors)	56.5 cm
Castor size	15 cm (6 inches)
Unit weight	6.7 kg (14.7 lbs)
Colour options	Blue, burgundy, purple, black, champagne
Warranty period	12 months Please note: Brake maintenance and adjustment are not covered under warranty