




**SKU:**

## **HERO MEDICAL 4- WHEELED MOBILITY WALKER - WALKING FRAME WITH ROLLATOR**

**\$169.99**

The Hero Medical Four Wheeled Seat Walker (Rollator) is much loved by ILS customers for every day use. Ideal for both indoor and outdoor use, the seat walker has a safe working limit of 140kg. Key features include a lightweight aluminium frame of 5kg, easy-to-operate functions and safety features such as a light touch lock brake system and large front wheels to minimise wheel wobble and increase grip. It easily folds without removing the storage bag and best of all, there are five colours to choose from.



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU    | Price    | Description | Colour    |
|----------------------------------------------------------------------------------|--------|----------|-------------|-----------|
|  | S11167 | \$169.99 |             | Burgundy  |
|  | S11991 | \$169.99 |             | Black     |
|  | S11990 | \$169.99 |             | Champagne |
|  | S11168 | \$169.99 |             | Blue      |
|  | S11169 | \$169.99 |             | Purple    |

## **PRODUCT DESCRIPTION**

The Hero Medical Four Wheeled Seat Walker (Rollator) is perfect for indoor and outdoor use for those with limited mobility. It is easy to transport and operate, and has many safety features.

Key features include:

- Lightweight aluminium frame of 5kg
- Folds easily
- Swivel front wheels for easy maneuvering
- 15cm (6 inches) flat ridged wheels to minimise wheel wobble and increase grip
- Adjustable handle bars
- Ergonomic hand grips
- Light touch hand lock brake system
- Padded easy fold rounded back rest
- Spill proof nylon padded seat
- Waterproof removable nylon bag with snaps

## ADDITIONAL INFORMATION

|                                             |                                                                                           |
|---------------------------------------------|-------------------------------------------------------------------------------------------|
| <b>Brand</b>                                | Hero Medical                                                                              |
| <b>Weight capacity / Max SWL</b>            | 140kg (308 lbs)                                                                           |
| <b>Overall height (floor to handles)</b>    | 77 cm – 92 cm                                                                             |
| <b>Overall length</b>                       | 69 cm (27 inches)                                                                         |
| <b>Overall width</b>                        | 62 cm (24.4 inches)                                                                       |
| <b>Seat height</b>                          | 55 cm (21.6 inches)                                                                       |
| <b>Seating width</b>                        | 44 cm (17.3 inches)                                                                       |
| <b>Seat depth</b>                           | 33 cm (12.9 inches)                                                                       |
| <b>Folded length</b>                        | 82cm (32.2 inches)                                                                        |
| <b>Folded width</b>                         | 23 cm (9 inches)                                                                          |
| <b>Outer width (including back castors)</b> | 56.5 cm                                                                                   |
| <b>Castor size</b>                          | 15 cm (6 inches)                                                                          |
| <b>Unit weight</b>                          | 6.7 kg (14.7 lbs)                                                                         |
| <b>Colour options</b>                       | Blue, burgundy, purple, black, champagne                                                  |
| <b>Warranty period</b>                      | 12 months<br>Please note: Brake maintenance and adjustment are not covered under warranty |