




SKU:

HERO MEDICAL 4- WHEELED MOBILITY WALKER - WALKING FRAME WITH ROLLATOR

\$169.99

The Hero Medical Four Wheeled Seat Walker (Rollator) is much loved by ILS customers for every day use. Ideal for both indoor and outdoor use, the seat walker has a safe working limit of 140kg. Key features include a lightweight aluminium frame of 5kg, easy-to-operate functions and safety features such as a light touch lock brake system and large front wheels to minimise wheel wobble and increase grip. It easily folds without removing the storage bag and best of all, there are five colours to choose from.



Scan QR code to
view product online

VARIATIONS





Image	SKU	Price	Description	Colour
	S11167	\$169.99		Burgundy
	S11991	\$169.99		Black
	S11990	\$169.99		Champagne
	S11168	\$169.99		Blue

Image	SKU	Price	Description	Colour
	S11169	\$169.99		Purple

PRODUCT DESCRIPTION

The Hero Medical Four Wheeled Seat Walker (Rollator) is perfect for indoor and outdoor use for those with limited mobility. It is easy to transport and operate, and has many safety features.

Key features include:

- Lightweight aluminium frame of 5kg
- Folds easily
- Swivel front wheels for easy maneuvering
- 15cm (6 inches) flat ridged wheels to minimise wheel wobble and increase grip
- Adjustable handle bars
- Ergonomic hand grips
- Light touch hand lock brake system
- Padded easy fold rounded back rest
- Spill proof nylon padded seat
- Waterproof removable nylon bag with snaps

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight capacity / Max SWL	140kg (308 lbs)
Overall height (floor to handles)	77 cm - 92 cm
Overall length	69 cm (27 inches)
Overall width	62 cm (24.4 inches)
Seat height	55 cm (21.6 inches)
Seating width	44 cm (17.3 inches)
Seat depth	33 cm (12.9 inches)
Folded length	82cm (32.2 inches)
Folded width	23 cm (9 inches)
Outer width (including back castors)	56.5 cm
Castor size	15 cm (6 inches)
Unit weight	6.7 kg (14.7 lbs)
Colour options	Blue, burgundy, purple, black, champagne
Warranty period	12 months Please note: Brake maintenance and adjustment are not covered under warranty