

×

SKU:

HERO MEDICAL 6" STANDARD SEAT WALKER

\$179.99

The Hero Medical 6" Standard Seat Walker is a compact, four-wheeled walking aid that offers safe, stable mobility support indoors and outdoors. With a padded seat, height-adjustable handles and built-in storage, it supports everyday independence.



Scan QR code to view product online



VARIATIONS

Image	SKU	Price	Description	Colours
×	S11167	\$179.99		Burgundy
×	S11991	\$179.99		Black
×	S11990	\$179.99		Champagne
×	S11168	\$179.99		Blue
×	S11169	\$179.99		Purple



PRODUCT DESCRIPTION

The Hero Medical 6" Standard Seat Walker is a practical mobility aid for individuals needing assistance with walking and resting. Its lightweight aluminium frame folds easily for transport and is fitted with 6" wheels to navigate a variety of surfaces, whether indoors or outdoors.

Height-adjustable handles help ensure ergonomic positioning, while light-touch loop brakes enhance control and stopping ease.

A cushioned seat and backrest support comfort during breaks, reducing fatigue and promoting safety. An under-seat storage bag provides convenient space for carrying personal belongings, making this walker a functional solution for daily use.

Who it's for

- For clients who need a cost-effective, standard walker to support mobility and reduce fatigue.
- For therapists prescribing a durable, compact walker with adjustable handbrakes to assist safe gait patterns.

Key features

- Lightweight aluminium frame easy to manoeuvre and store
- Four-wheeled design offers stability and smooth movement
- Padded seat and backrest support rest and comfort during use
- Adjustable handle height promotes correct posture and safe walking
- Light-touch loop brakes improve safety and control
- Folding frame with storage bag allows easy transport and keeps items accessible



ADDITIONAL INFORMATION

Safe working load (SWL) 140 kg

Gross Weight 6.7 kg

Seat Height 55 cm

Seat Width 44 cm

Handle Height 77 cm - 92 cm

Outer Width 56.5 cm

Front Castor 6"

Rear Castor 6"

Frame Material Aluminium

Brand Hero Medical