



**SKU:**

## **HERO MEDICAL 6" STANDARD SEAT WALKER**

**\$179.99**

The Hero Medical 6" Standard Seat Walker is a compact, four-wheeled walking aid that offers safe, stable mobility support indoors and outdoors. With a padded seat, height-adjustable handles and built-in storage, it supports everyday independence.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S11167	\$179.99		Burgundy
	S11991	\$179.99		Black
	S11990	\$179.99		Champagne
	S11168	\$179.99		Blue
	S11169	\$179.99		Purple

## PRODUCT DESCRIPTION

The Hero Medical 6" Standard Seat Walker is a practical mobility aid for individuals needing assistance with walking and resting. Its lightweight aluminium frame folds easily for transport and is fitted with 6" wheels to navigate a variety of surfaces, whether indoors or outdoors.

Height-adjustable handles help ensure ergonomic positioning, while light-touch loop brakes enhance control and stopping ease.

A cushioned seat and backrest support comfort during breaks, reducing fatigue and promoting safety.

An under-seat storage bag provides convenient space for carrying personal belongings, making this walker a functional solution for daily use.

### Who it's for

- For clients who need a cost-effective, standard walker to support mobility and reduce fatigue.
- For therapists prescribing a durable, compact walker with adjustable handbrakes to assist safe gait patterns.

### Key features

- **Lightweight aluminium frame** – easy to manoeuvre and store
- **Four-wheeled design** – offers stability and smooth movement
- **Padded seat and backrest** – support rest and comfort during use
- **Adjustable handle height** – promotes correct posture and safe walking
- **Light-touch loop brakes** – improve safety and control
- **Folding frame with storage bag** – allows easy transport and keeps items accessible

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Weight capacity / Max SWL</b>	140kg (308 lbs)
<b>Overall height (floor to handles)</b>	77 cm – 92 cm
<b>Overall length</b>	69 cm (27 inches)
<b>Overall width</b>	62 cm (24.4 inches)
<b>Seat height</b>	55 cm (21.6 inches)
<b>Seating width</b>	44 cm (17.3 inches)
<b>Seat depth</b>	33 cm (12.9 inches)
<b>Folded length</b>	82cm (32.2 inches)
<b>Folded width</b>	23 cm (9 inches)
<b>Outer width (including back castors)</b>	56.5 cm
<b>Castor size</b>	15 cm (6 inches)
<b>Unit weight</b>	6.7 kg (14.7 lbs)
<b>Colour options</b>	Blue, burgundy, purple, black, champagne
<b>Warranty period</b>	12 months Please note: Brake maintenance and adjustment are not covered under warranty