



SKU:

HERO MEDICAL TRAVEL LITE PORTABLE SEAT WALKER

\$349.99

The Hero Medical Travel Lite Portable Seat Walker provides lightweight, foldable support for daily mobility. With 8" wheels, ergonomic handles and seated comfort, it's suited for both indoor and outdoor use.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Colours
S11204	\$349.99		Blue
S12385	\$349.99		Silver

GALLERY IMAGES



PRODUCT DESCRIPTION

The Hero Medical Travel Lite Portable Seat Walker is a compact mobility aid designed to support independence and mobility in community settings. Its lightweight aluminium frame folds easily, making it ideal for transport and storage in vehicles or smaller spaces.

Equipped with 8" wheels, it glides smoothly over indoor and outdoor surfaces, while a padded seat and backrest offer support during rest breaks.

Height-adjustable ergonomic handles encourage optimal posture, and light-touch loop brakes provide controlled stopping and secure parking.

Additional features include a removable storage bag and a crutch holder, enhancing day-to-day convenience and personal independence.

Who it's for

- For clients seeking a lightweight, portable seat walker for everyday mobility.
- For therapists prescribing a compact, travel-ready walker with ergonomic features and durable construction.

Key features

- **Folding aluminium frame** – supports portability and easy storage
- **8" wheels** – enable smooth movement across indoor and outdoor surfaces
- **Padded seat and backrest** – provide comfortable rest during outings
- **Height-adjustable ergonomic handles** – support safe and functional walking posture
- **Light-touch loop brakes** – ensure user-controlled stopping and secure parking
- **Removable storage bag and crutch holder** – promote convenience and independence
- **Compact folding mechanism** – fits into car boots and small storage areas

ADDITIONAL INFORMATION

Safe working load (SWL)	140 kg
Gross Weight	8 kg
Seat Height	55 cm
Seating Width	46 cm
Handle Height	78 cm – 94 cm
Outer Width	60 cm
Front Castor	8"
Rear Castor	8"
Frame Material	Aluminium
Colours	Blue, Silver
Brand	Hero Medical

[Hero Medical Outdoor Lite Portable Seat Walker User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

HOW DO I ADJUST THE HEIGHT?

For the Travel Lite Portable Outdoor Seat Walker, the height is adjustable on the two side bars. There are height adjustment locking pins. You will need to push in the pin to move it on to the next height level.

Overall Height

Min 78cm
Max 90cm

WHAT IS THE FOLDED DIMENSIONS?

Folded width is 27cm.