



SKU:

## **TRAVEL LITE PORTABLE OUTDOOR SEAT WALKER WITH SEAT AND BAG + CRUTCH HOLDER**

**\$349.99**

The Travel Lite Portable Outdoor Seat Walker is a new modern seat walker design by Hero. It is designed to give maximum mobility for the user with many important features such as the quick adjust handle height and easy fold system, superior manoeuvrability and seating comfort. It folds into a compact unit for ease of transportation and is low maintenance. International Standards tested.



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Colours
 A blue and black HERO wheelchair with a large black storage bag attached to the front. The brand name 'HERO' is visible on the backrest.	S11204	\$349.99		Blue
 A silver and black HERO wheelchair with a large black storage bag attached to the front. The brand name 'HERO' is visible on the backrest.	S12385	\$349.99		Silver

## **PRODUCT DESCRIPTION**

### **Travel Lite Portable Outdoor Seat Walker**

This Travel Lite Portable Outdoor Seat Walker with Seat and Bag + Crutch Holder by Freedom Healthcare offers the following features:

- Woven padded nylon sling Backrest
- Light-touch lock brakes, ideal for arthritis sufferer
- Woven breathable nylon seat
- Dual purpose removable shoulder bag
- Easy fold front 8 inch Large outdoor castors wheels
- Foldable rear legs with knob adjustment system
- Handle bars with directional adjustable for height with push button adjustment system
- Walking stick attachment included

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Weight capacity / Max SWL</b>	140kg
<b>Overall height (floor to handles)</b>	Min 78cm Max 94cm
<b>Overall length</b>	69cm
<b>Overall width</b>	64cm
<b>Seat depth</b>	25cm
<b>Seating width</b>	46cm
<b>Seat height</b>	55cm
<b>Folded width</b>	27cm
<b>Outer width (including back castors)</b>	60cm
<b>Castor size</b>	20cm (8 inches)
<b>Unit weight</b>	8kg
<b>Colours</b>	Blue, Silver
<b>Warranty</b>	12 Months Warranty (Brake maintenance and adjustment is not covered under warranty.)