



SKU:

HERO MEDICAL TRAVEL LITE PORTABLE SEAT WALKER

\$349.99

The Hero Medical Travel Lite Portable Seat Walker provides lightweight, foldable support for daily mobility. With 8" wheels, ergonomic handles and seated comfort, it's suited for both indoor and outdoor use.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S11204	\$349.99		Blue
	S12385	\$349.99		Silver

PRODUCT DESCRIPTION

The Hero Medical Travel Lite Portable Seat Walker is a compact mobility aid designed to support independence and mobility in community settings. Its lightweight aluminium frame folds easily, making it ideal for transport and storage in vehicles or smaller spaces.

Equipped with 8" wheels, it glides smoothly over indoor and outdoor surfaces, while a padded seat and backrest offer support during rest breaks.

Height-adjustable ergonomic handles encourage optimal posture, and light-touch loop brakes provide controlled stopping and secure parking.

Additional features include a removable storage bag and a crutch holder, enhancing day-to-day convenience and personal independence.

Who it's for

- For clients seeking a lightweight, portable seat walker for everyday mobility.
- For therapists prescribing a compact, travel-ready walker with ergonomic features and durable construction.

Key features

- **Folding aluminium frame** – supports portability and easy storage
- **8" wheels** – enable smooth movement across indoor and outdoor surfaces
- **Padded seat and backrest** – provide comfortable rest during outings
- **Height-adjustable ergonomic handles** – support safe and functional walking posture
- **Light-touch loop brakes** – ensure user-controlled stopping and secure parking
- **Removable storage bag and crutch holder** – promote convenience and independence
- **Compact folding mechanism** – fits into car boots and small storage areas

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight capacity / Max SWL	140kg
Overall height (floor to handles)	Min 78cm Max 94cm
Overall length	69cm
Overall width	64cm
Seat depth	25cm
Seating width	46cm
Seat height	55cm
Folded width	27cm
Outer width (including back castors)	60cm
Castor size	20cm (8 inches)
Unit weight	8kg
Colours	Blue, Silver
Warranty	12 Months Warranty (Brake maintenance and adjustment is not covered under warranty.)