



SKU:

TRAVEL LITE PORTABLE OUTDOOR SEAT WALKER WITH SEAT AND BAG + CRUTCH HOLDER

\$349.99

The Travel Lite Portable Outdoor Seat Walker is a new modern seat walker design by Hero Medical. It is designed to give maximum mobility for the user with many important features such as the quick adjust handle height and easy fold system, superior manoeuvrability and seating comfort. It folds into a compact unit for ease of transportation and is low maintenance. International Standards tested.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S11204	\$349.99		Blue
	S12385	\$349.99		Silver

PRODUCT DESCRIPTION

Travel Lite Portable Outdoor Seat Walker

This Travel Lite Portable Outdoor Seat Walker with Seat and Bag + Crutch Holder by Freedom Healthcare offers the following features:

- Woven padded nylon sling Backrest
- Light-touch lock brakes, ideal for arthritis sufferer
- Woven breathable nylon seat
- Dual purpose removable shoulder bag
- Easy fold front 8 inch Large outdoor castors wheels
- Foldable rear legs with knob adjustment system
- Handle bars with directional adjustable for height with push button adjustment system
- Walking stick attachment included

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight capacity / Max SWL	140kg
Overall height (floor to handles)	Min 78cm Max 94cm
Overall length	69cm
Overall width	64cm
Seat depth	25cm
Seating width	46cm
Seat height	55cm
Folded width	27cm
Outer width (including back castors)	60cm
Castor size	20cm (8 inches)
Unit weight	8kg
Colours	Blue, Silver
Warranty	12 Months Warranty (Brake maintenance and adjustment is not covered under warranty.)