

SKU: S13001



DAYS 2-IN-1 SEAT WALKER/WHEELCHAIR

\$299.99

The Days 2-in-1 Seat Walker/Wheelchair offers versatile support for those with limited or fluctuating mobility. It functions as both a walking aid and a lightweight wheelchair, providing freedom to walk or be assisted when fatigued.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Days 2-in-1 Seat Walker/Wheelchair is a practical mobility solution for individuals who can walk short distances but may need seated support during longer outings.

It enables greater independence by allowing supported walking but also transitions to seated wheelchair mobility when required.

The frame converts easily from seat walker to transit wheelchair. It features a cushioned seat and canvas backrest for comfort, ergonomic handgrips with loop brakes for safe handling, and swing-away footrests for stable positioning during use as a wheelchair.

The lightweight aluminium frame folds for transport and storage, while large 7.5" castors support use across indoor and outdoor surfaces.

Who it's for

- For clients who can walk short distances with support but benefit from a seated mobility option when tired.
- For therapists supporting clients with variable endurance who need both walking support and assisted transport.

Key features

- **Dual-use design** – functions as both a seat walker and transit wheelchair
- **Foldable aluminium frame** – strong, lightweight and easy to store or transport
 - **Padded seat and canvas backrest** – for comfort when seated
 - **Ergonomic loop brakes** – allow safe control while walking
- **Swing-away footrests** – support comfort and safety in wheelchair mode
- **7.5” castors** – suitable for smooth travel on indoor and outdoor surfaces

ADDITIONAL INFORMATION

Safe Working Load (SWL)	160 kg
Gross Weight	8.8 kg
Seat Height	56 cm
Seat Width	35 cm
Handle Height	80 cm – 93 cm
Outer Width	58.5 cm
Front Castor	7.5"
Rear Castor	7.5"
Frame Material	Aluminium
Brand	Performance Health

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

WHAT IS THE SIDE-TO-SIDE WIDTH OF THE SEAT?

For the Seat Walker/Wheelchair - Days Transit 2 in 1, the seat width is 350mm across

WHAT IS YOUR RETURNS POLICY?

Please refer to return policy <https://ilsau.com.au/return-policy/>