

SKU: S13001




DAYS 2 IN 1 TRANSIT ROLLATOR - WALKER WITH SEAT

\$299.99

The Days 2 in 1 Transit Rollator combines the features of a Walker or Rollator and a Transit / Transport Chair. Users are able to operate it as a superior four-wheeled rollator on its own or even just simply flip down the footrest and convert it into a transit chair when needed. The Days 2 in 1 rollator provides users and carers a choice on how to settle down on the ultimate comfort of the patient.

We recommend [Wheelchair/Rollator - FUSION 2 IN 1 - Hero](#)



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Days 2 in 1 transit rollator combines the features of a walker/rollator and a transit /transport chair into one unit. You can use it as a deluxe four wheeled rollator by itself or when you are tired, you can just simply engage the flip down footrest and convert into a transit chair and allow your carer to push you around.

It has large 7.5" caster wheels for a smooth ride, perfect for indoor and outdoor use. Handles adjust to the user's height while the reliable, loop style brakes control speed and lock wheels when seated.

Key features include:

- Large 7.5" wheels with wheel locks
- Fold down footrests make transporting a patient easy
- Locking, easy-to-use loop style brakes
- Padded seat with easy-folding system
- Weight capacity: 160kg
- Available in red only.

ADDITIONAL INFORMATION

Brand	Days Healthcare
Max SWL	160kg (352.74 lbs)
Seat height	56 cm (22 inches)
Seat width	35 cm (13.7 inches)
Seat depth	33 cm (12.9 inches)
Overall height	80 - 93 cm (31.4-36.6 inches)
Overall width	58.5 cm (23 inches)
Overall depth	74 cm (29.1 inches)
Weight	8.8 kg (19.4 lbs)
Wheel	7.5"
Color	Red
Warranty	12 months