



**SKU:** S11201

## **SELF HELP POLE - FREE STANDING**

**\$299.99**

The Royale Self Help Pole is designed to provide a practical and supportive solution to assist individuals with mobility challenges. Makes it easier for users to reposition themselves, sit up and transfer in and out of bed, via a secure and stable support handle.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

This freestanding self help pole helps promote independence by allowing individuals to reposition themselves in bed via a convenient overhead support handle. The support handle has an adjustable strap to either lower or increase the height. The self-help pole can be positioned on any type of surface and is fitted suction rubber tips.

- Dismantles into two pieces
- Sits at the bed head or can be positioned beside the bed
- Support handle has a height adjustment strap
- Easy to install

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Overall Height</b>	183 cm
<b>Overall Length</b>	80 cm
<b>Overall Width</b>	81 cm
<b>Unit Weight</b>	11.3 kg
<b>Max SWL</b>	140 kg