




SKU: S087

INVACARE PORTABLE FOLDING WHEELCHAIR RAMP - SENIOR

~~\$2,999.99~~ \$2,999.99

A simple concertina folding action enables the Invacare Portable Ramp to be deployed in seconds. As with the personal and multi-purpose models each is remarkably lightweight and strong in use. The Senior Range offers you a flexible variety of ramps to cope with obstacle. Passes Australian Standards AS3856-1.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Length
	S087	\$699.99		90 cm
	S1660	\$799.99		120 cm
	S081	\$899.99		135cm
	S090	\$949.99		150 cm
	S085	\$999.99		165 cm
	S7645	\$1,199.99		185 cm
	S082	\$1,499.99		200 cm
	S7646	\$1,999.99		250 cm
	S7647	\$2,499.99		300 cm
	S1691	\$2,999.99		350 cm

PRODUCT DESCRIPTION

- **Portable** - Invacare Portable Ramps can be quickly folded into a single compact unit - allowing for easy storage or transport in any vehicle or home environment.
- **Lightweight** - Made from a unique combination of unidirectional glass and carbon fibre materials. Invacare Portable Ramps are very strong, yet light enough to be portable.
- **Compact** - Multi-paneled design held together with a patented hinge system connected by webbing enabling Invacare Portable Ramps to fold with a concertina action, allowing for a much smaller width when folded.
- **Safe** - Clearly marked and grip tape covered, Invacare Portable Ramps are also fitted with pads at the top and base. Invacare Portable Ramps provide a non-slip surface, for all weather conditions.

Ramp Sizes

Length	Width	Weight
90cm	75cm	4.8kg
120cm	75cm	8kg
135cm	75cm	NA
150cm	75cm	NA
165cm	75cm	10.8kg
185cm	75cm	NA
200cm	75cm	11.7kg
250cm	75cm	NA
300cm	75cm	21.1kg
350cm	75cm	24.7kg

ADDITIONAL INFORMATION

Brand	Invacare
Max height capacity	22.5 cm
Ramp width	75 cm
Safe working load	300kg