


SKU:

INVACARE DECPAC PORTABLE SENIOR RAMP

~~\$3,699.99~~ **\$3,699.99**

The Invacare Decpac Portable Senior Ramp is a wide, foldable ramp designed for secure access across steps and thresholds. Its larger size supports broader mobility devices while remaining lightweight and portable.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Select size
	S087	\$869.99		90 cm
	S1660	\$969.99		120 cm
	S081	\$999.99		135 cm
	S090	\$1,069.99		150 cm

Image	SKU	Price	Description	Select size
	S085	\$1,169.99		165 cm
	S7645	\$1,599.99		185 cm
	S082	\$1,899.99		200 cm
	S7646	\$2,599.99		250 cm

Image	SKU	Price	Description	Select size
	S7647	\$2,999.99		300 cm
	S1691	\$3,699.99		350 cm

PRODUCT DESCRIPTION

This Invacare Decpac Portable Senior Ramp features a concertina folding mechanism, allowing it to be deployed quickly without assembly. Its multi-panel design, held together with a patented hinge system connected by webbing, enables compact folding for easy storage and transport. Constructed from a combination of uni-directional glass and carbon fibre materials, the ramp offers strength while remaining lightweight. The surface is covered with grip tape and fitted with pads at the top and base to provide a non-slip surface in all weather conditions.

Who it's for

- For clients using wheelchairs or powered devices who require assistance navigating steps or raised surfaces.
- For therapists seeking portable ramp solutions for clients to access various environments safely.

Key features

- **Compact foldable design** – convenient for storage and transport
- **Lightweight carbon and fibreglass build** – strong yet manageable
- **Concertina folding mechanism** – simplifies storage and transport
- **Grip tape surface and padded ends** – promotes safety on varied terrain
- **Certified to AS/NZ standards** – supports up to 300 kg
- **No setup required** – folds out instantly for use

ADDITIONAL INFORMATION

Brand	Invacare
Max height capacity	22.5 cm
Ramp width	75 cm
Safe working load	300 kg