



SKU: S2445

SHOWER CHAIR - HEAVY DUTY

\$399.99

Height adjustable shower chair with arms and removable backrest. Non-padded seat and backrest. Comfortable and lightweight for easy transportation and manoeuvring. Wide and comfortable for both standard and bariatric market. Aluminium and durable frame.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Key Features

- Height adjustable shower chair with removable arms and backrest.
- Comfortable and lightweight for easy transportation and storage
- Wide and comfortable for both standard and bariatric clients
- Rubber non-slip tips on feet for added safety
- Lightweight aluminium frame

ADDITIONAL INFORMATION

| | |
|---|-------------------|
| Brand | Safety & Mobility |
| Overall Width | 67 cm |
| Overall Depth | 43 cm |
| Seat Width | 52 cm |
| Width Between the Arms | 62 cm |
| Seat Height Range | 48 cm – 58 cm |
| Armrest Height (from seat to top of arm) | 15 cm |
| Weight Capacity / Max SWL | 150 kg |