




**SKU:** S11195

## **HERO MEDICAL STANDARD SHOWER CHAIR**

**\$149.99**

The Hero Medical Standard Shower Chair provides reliable, supportive seating for people who need seated stability while showering. With its rust-resistant aluminium frame, contoured seat and fixed armrests, it's a practical choice for home and aged care settings.



Scan QR code to  
view product online

### **GALLERY IMAGES**



### **PRODUCT DESCRIPTION**

The Hero Medical Standard Shower Chair supports safe, comfortable seated bathing for individuals with limited strength, balance or mobility. The lightweight aluminium frame is durable and rust-resistant, with a safe working load of 205 kg.

Height-adjustable legs allow the seat height to be customised to suit different individuals and shower setups. The contoured seat features drainage holes to promote hygiene and reduce water pooling.

Fixed raised armrests assist with safe sit-to-stand transfers, providing extra security during use.

Non-slip rubber feet help keep the chair stable on wet shower floors.

### Who it's for

- For clients who need seated support during daily showering due to risk of falls.
- For therapists recommending reliable shower seating for clients in home or aged care settings.

### Key features

- **45 cm aluminium frame** – lightweight, rust-resistant and durable
- **Height-adjustable legs** – customise seat height to suit individual needs
- **Contoured seat with drainage holes** – enhances hygiene and comfort
  - **Fixed raised armrests** – assist with safe sit-to-stand transfers
  - **Non-slip rubber feet** – improve grip and stability on wet floors
- **Safe working load of 205 kg** – tested and compliant with Australian Standards

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	205 kg
<b>Gross Weight</b>	3 kg
<b>Overall Width</b>	54 cm
<b>Overall Depth</b>	54 cm
<b>Seat Height</b>	42 cm – 54 cm
<b>Seat Width</b>	44 cm
<b>Width Between the Arms</b>	45 cm
<b>Armrest Height</b>	22 cm
<b>Brand</b>	Hero Medical

[Hero Medical Standard Shower Chair User Manual](#)

### **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)