




**SKU:** S11195

## **SHOWER CHAIR - ALUMINIUM RUST FREE**

**\$149.99**

This Hero Medical shower chair is made from a lightweight aluminium frame (3kg) and has a comfortable contoured plastic seat. The seat has drainage holes and provides posture support whilst showering. This shower chair has a safe usage weight of 205kg which has been tested to Australian standards.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

Shower chairs are an ideal safety measure for users with mobility limitations as it allows them to be seated whilst showering and reduces the risk of falling. The shower chair comes with height adjustable legs, allowing multiple seat height options that are suitable for every user. The shower chair has a contoured plastic seat with drainage holes, anti-slip tips to ensure stability and the frame has been treated with corrosion resistance and powder coated paint.

### **Key Features**

- Lightweight aluminium frame
- Height adjustable legs
- Raised arm supports to allow easier transfer on/off
- Durable anti-slip tips
- SWL 205kg
- 45cm width between the arms

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Frame Material</b>	Aluminium
<b>Overall Width</b>	54cm
<b>Overall Depth</b>	54cm
<b>Seat Width</b>	45cm
<b>Width Between the Arms</b>	46cm
<b>Seat Height Range</b>	42cm – 54cm
<b>Armrest Height (from seat to top of arm)</b>	22cm
<b>Unit Weight</b>	3kg
<b>Weight Capacity / Max SWL</b>	205kg