



SKU: S11195

HERO MEDICAL STANDARD SHOWER CHAIR

\$149.99

The Hero Medical Standard Shower Chair provides reliable, supportive seating for people who need seated stability while showering. With its rust-resistant aluminium frame, contoured seat and fixed armrests, it's a practical choice for home and aged care settings.



Scan QR code to view product online



PRODUCT DESCRIPTION

The Hero Medical Standard Shower Chair supports safe, comfortable seated bathing for individuals with limited strength, balance or mobility. The lightweight aluminium frame is durable and rust-resistant, with a safe working load of 205 kg.

Height-adjustable legs allow the seat height to be customised to suit different individuals and shower setups. The contoured seat features drainage holes to promote hygiene and reduce water pooling. Fixed raised armrests assist with safe sit-to-stand transfers, providing extra security during use. Non-slip rubber feet help keep the chair stable on wet shower floors.

Who it's for

- For clients who need seated support during daily showering due to risk of falls.
- For therapists recommending reliable shower seating for clients in home or aged care settings.

Key features

- 45 cm aluminium frame lightweight, rust-resistant and durable
- Height-adjustable legs customise seat height to suit individual needs
- Contoured seat with drainage holes enhances hygiene and comfort
- Fixed raised armrests assist with safe sit-to-stand transfers
- Non-slip rubber feet improve grip and stability on wet floors
- Safe working load of 205 kg tested and compliant with Australian Standards



ADDITIONAL INFORMATION

Safe Working Load (SWL) 205 kg

Gross Weight 3 kg

Overall Width 54 cm

Overall Depth 54 cm

Seat Height 42 cm - 54 cm

Seat Width 44 cm

Width Between the Arms 45 cm

Armrest Height 22 cm

Brand Hero Medical