




SKU: S2229

SHOWER COMMODE - SELF PROPELLED FOLDING MOBILE 445MM WITH SWING AWAY LEGRESTS

From: \$2,499.99

The Self Propelled Folding Mobile Shower Commode assists in the transportation and care for patient showering and toileting. This commode comes with a 3-strap backrest for comfort and easy access. Also comes with an integrated push-handle for manoeuvrability and has swing-up arms for easier side transfers. Equipped with swing-away, height adjustable footrests, 2x125 mm total locking castors and 2 x 22" self-propelled wheels with push-to-lock brakes. This commode has been independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.

Please note, seat sold separately



Scan QR code to
view product online

PRODUCT DESCRIPTION

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- 3 strap back rest for comfort and easy access
- Integrated push-handle for manoeuvrability
- Swing-up arms for unrestricted side transfer
- Swingaway and removable height adjustable footrest
- 2×125 mm total locking castors
- 2×22"self-propel wheels inc. push to lock brakes
- Push rims available as optional extra
- Padded seat not included

ADDITIONAL INFORMATION

Overall Width	630 mm
Overall Depth	1050 mm
Overall Height	1005 mm
Seat Width	445 mm
Seat Depth	430 mm
Floor To Seat Height	550 mm
Armrest Type	Loop
Back Height	400 mm
Footrest Type	Swing Away
Product Weight	15 kg
Max User Weight	200 kg
Brand	KCare