




SKU: S2229

SHOWER COMMODE - MOBILE - SELF PROPELLED - FOLDING

From: \$2,499.99

The folding shower commode assists with the transportation of a patient for showering and toileting. It is easily positioned over the toilet for toileting and folds neatly for transportation or storage where space is limited. The padded commode seat (available as open or closed front) features adjustable height swing-away leg-rests assist patient transfers as do it's swing-back arms for side transfer. The hinge point is in the middle of the seat base of the shower commode allowing it to fold inwards. Australian made.

Seat sold separately



Scan QR code to
view product online

PRODUCT DESCRIPTION

A removable PVC slatted backrest aids in storage and patient access and the folding shower commode moves on a pair 2 x 125 mm castors and 2 x 22" rear wheels. Built with a durable stainless steel frame. Padded Seat is optional.

ADDITIONAL INFORMATION

Seat Width	445 mm
Seat Depth	430 mm
Seat Height	590 mm
Arm Height	240 mm
Back Height	365 mm
Seat to Footplate	380 mm – 520 mm
Overall Width	610 mm
Max SWL	120 kg
daily-living	Commodes