


**SKU:** S11968



## **HERO MEDICAL ALUMINIUM RUST FREE PADDED SHOWER STOOL AND CHAIR**

**\$169.99**

This is the Hero Medical padded shower stool which is a great option for supporting individuals experiencing difficulties when showering and washing. This product has a safe usage load of 205kg and has been tested and is compliant with Australian Standards. The product has a lightweight sturdy aluminium frame which has been treated with corrosive resistance powder paint. The seat is padded and provides comfort for users. The frame and seat are easily cleaned and maintained.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

Shower stools are an excellent falls prevention measure for individuals with limited mobility to improve safety. The Hero shower stool can assist individuals maintaining independence as it enables users to be seated whilst showering. The product comes with height adjustable legs which will alter the seat height to the users needs, raised armrests for additional stability when raising and lowering and a padded seat. The shower chair is equipped with anti-slip rubber tips and has been treated with corrosive resistance powder paint. This model has a SWL of 205kg and is 46cm in width between the arms.

### **Key Features**

- Lightweight aluminium frame
- Padded seat
- Height adjustable legs
- Raised armrests for easy transferring on/off
- Durable anti-slip tips
- Frame has been treated with corrosive resistance powder paint
- SWL 205kg

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Frame</b>	Aluminium
<b>Overall Width</b>	60cm
<b>Overall Depth</b>	47cm
<b>Seat Width</b>	35cm
<b>Width Between the Arms</b>	46cm
<b>Seat Height Range</b>	46cm -59cm
<b>Armrest Height (from seat to top of arm)</b>	22cm
<b>Unit Weight</b>	4kg
<b>Weight Capacity / Max SWL</b>	205kg