

SKU:

SLIDE SHEET PREMIUM

~~\$69.99~~ **\$69.99**

The Slide Sheet Premium is silicon coated on both sides and is very easy to use. It slips consistently in all directions and is made of soft polyester taffeta. The Slide Sheet will retain it's slip for over 150 washes and is machine washable up to 75c.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S1529	\$59.99		2.0m x 1.5m
	S10348	\$69.99		3.0m x 1.5m
	S10347	\$39.99		1.0m x 1.5m

PRODUCT DESCRIPTION

You can tumble dry at cool temperature or hang to dry for longer life. SWL of 360 kg (Australian Wool Testing Authority). TGA Registered.

ADDITIONAL INFORMATION

Brand

Haines Medical

Sizes

3.0m x 1.5m, 2.0m x 1.5m, 1.0m x 1.5m