



SKU:

SLIDE SHEET PREMIUM

~~\$69.99~~ **\$69.99**

These premium slide sheets are washable for reuse with multiple patients. Soft and quiet, they slip extremely well in all directions and are the slide sheets of choice for many of Australia's leading manual handling trainers. Made from soft, high-grade silicone coated polyester taffeta (slippery on both sides). Are washable, reusable, radiolucent and have a SWL of up to 360 kg.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S10347	\$39.99		1.0m x 1.45m
	S1529	\$59.99		2.0m x 1.45m
	S10348	\$69.99		3.0m x 1.5m

PRODUCT DESCRIPTION

- Soft, high-grade, silicone-coated polyester taffeta – Slippery on both sides (no wrong way up)
- Sheets are proven to retain their slip for up to 250 washes at 75°C (tumble dry on cool setting or hang to dry for longer life)
- Radiolucent, can be used for X-rays/CT scans
- Safe Working Load (SWL) up to 360 kg
- Reusable/Washable

ADDITIONAL INFORMATION

Brand	SlipperySally
Sizes	1.0m x 1.45m, 2.0m x 1.45m, 3.0m x 1.5m
Device Type	Small Assistive Patient Transfer Aids