




SKU:

SOLARI MEDICAL COCCYX WEDGE CUSHION

\$69.99

The Solari Medical Coccyx Wedge Cushion is designed to relieve pressure on the coccyx and lower spine, supporting improved posture and comfort for extended sitting periods. Its wedge shape and high-resilience foam construction provide consistent support for everyday use.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	SIZES
	S10353	\$69.99	Cushion - Coccyx Wedge Large Solari	Large 16" x 16"
	S10352	\$69.99	Cushion - Coccyx Wedge Standard Solari	Standard 16" x 14"

PRODUCT DESCRIPTION

The Solari Medical Coccyx Wedge Cushion features a wedge-shaped design that tilts the pelvis forward, promoting optimal spinal alignment and reducing strain on the lower back. The cutaway tailbone section helps alleviate pressure on the coccyx, offering targeted relief for those managing tailbone discomfort or lower back pain.

Constructed from high-resilience foam, this cushion delivers consistent support and durability, making it a reliable choice for long-term use. Lightweight and portable, it's suitable for a range of seating surfaces, including wheelchairs, office chairs and car seats, ensuring comfort and support wherever it's needed.

Who it's for

- For clients experiencing coccyx or lower back pain, or recovering from tailbone injuries.
- For therapists seeking a supportive, easy-to-maintain cushion that promotes pressure relief and optimal sitting posture.

Key features

- **Wedge-shaped design** – promotes forward pelvic tilt for better spinal alignment
- **Cutaway tailbone section** – reduces pressure on the coccyx, alleviating discomfort
- **High-resilience foam construction** – provides consistent support and durability
- **Lightweight and portable** – suitable for various seating environments

ADDITIONAL INFORMATION

SIZES	Standard 16" x 14", Large 16" x 16"
Low Height	1 cm
High Height	7 cm, 9 cm
Brand	Solari Medical
SWL	Unlimited
Risk Rating Filter	Low