




ABLE LIFE SPACE SAVER ROLLATOR

~~\$549.99~~ \$299.99

The World's Most Compact and Portable Rollator




Able life's new space saver rollator is the perfect unit to keep you independent and mobile on a daily basis. This stylish and unique rollator is a patented design that allows you to fold up to a 13-inch diameter. The rollator is ideal for traveling, grocery shopping, dining out or even just maneuvering along the small spaces around your home.



SKU: N/A

Category: [Walking Aids](#)

VARIATIONS

Image	SKU	Price	Description	Colour
	S17148	\$549.99 <u>\$299.99</u>		Walnut Black
	S17147	\$549.99 <u>\$299.99</u>		Blue
	S17149	\$549.99 <u>\$299.99</u>		Regal Rose

GALLERY IMAGES





PRODUCT DESCRIPTION

Portable and Lightweight

The rollator's unique folding abilities and its lightweight frame (at just 6.3kg) brands the Space Saver Rollator as one of the most portable walkers across the world. It can even be folded into the front seat of a car!

Strong and Durable

Regardless of its lightweight frame, the rollator can carry up to 115kg. The Space Saver Rollator is created with high-grade aircraft aluminum and is the perfect blend of durability and lightness.

The Space Saver Rollator is equipped with a full securing brake system with revisable brake pads

Portable – Swiftly and easily folds up for storage in smaller spaces

Convenient seat – Perfect for relaxation

Adjustable height – And carries up to 115kg

Modern look – Three trendy colors to choose from

6" wheels – Allows for easy maneuvering over any surface

ADDITIONAL INFORMATION

Colour	Walnut Black, Blue, Regal Rose
Weight	6.3kg
Max SWL	115kg
Collapsed Dimension	35.75"(91cm) X 13" (33cm) x 10" (25.4cm)
Height Adjustment	33"(84cm) - 36"(91cm)
Open Depth	23" (58cm)
Open Width	27.5" (70cm)
Open Handle Width	17" (43cm)
Seat Dimension	18.5 (47cm)" W X 5" (13cm)D
Brand	Able Life