




SKU:

SPECIAL TOMATO EIO PUSH CHAIR

The EIO Push Chair by Special Tomato is designed for children from early intervention to pre-adolescence and provides mild to moderate positioning support, with removable hip and chest supports included as standard.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The 3 point hip belt and chest harness with chest strap have a comfort cover while providing secure positioning. It's ready to go with sun canopy, underseat basket also standard accessories. For children who need to sleep while on the go, the backrest reclines 75 degrees. The stroller is easy to fold and features a lightweight aluminium frame. A raincover is also available. The size 1 and 2 Special Tomato Sitters and Seat Liners can be also be used in the stroller to provide extra postural support.

- Available in one size
- Max weight 40 kg

Accessories

- Rain cover

