


**SKU:** S12650

## STAINLESS STEEL PLATE SURROUND

**\$29.99**



The plate surround is an adjustable guard that can be fitted to any plate to assist individuals with one handed meal consumption. It is designed with clips attached to the rim that produce an inward sloping movement to assist with collection of food, and also prevent spillage.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

### Key Features

- Adjustable plate surround suitable for most types of plate
- Assists user for one handed meal consumption
- Dishwasher safe

## ADDITIONAL INFORMATION

<b>Brand</b>	ILS ADL
<b>Material</b>	Stainless steel
<b>Weight</b>	105gram