

SKU: S2296

## FORTE STANDARD HOSPITAL MATTRESS

~~\$699.99~~ **\$699.99**

The Forte Standard Hospital Mattress is ideal for low care applications where support, comfort, hygiene, durability, and stability is required for clients who are at low risk of pressure injury. This mattress comes with a Advantiflex waterproof cover which is a high-performance, 4-way stretch polyurethane material.



Scan QR code to  
view product online



## VARIATIONS

Image	SKU	Price	Description	Sizes
	S2296	\$499.99		Single
	S12254	\$699.99		King Single

## PRODUCT DESCRIPTION

The Forte Standard Hospital Mattress is a versatile, low-risk pressure care mattress designed for individuals requiring comfort and support in various care settings. Its durable construction and Advantiflex cover make it suitable for homecare, aged care and short-term hire applications. The mattress features a supportive foam core which provides comfort and support for repositioning and basic patient care. The mattress is fitted with an Advantiflex PU cover which is a 4-way stretch material for shear protection, immersion and envelopment. The cover is vapour permeable and breathable, resistant to bodily fluids, flexible and gentle to the skin and meets all related compliance requirements.

### Who it's for

- For clients with low pressure care risk needing a hygienic, supportive mattress.
- For therapists seeking an affordable, durable and supportive static foam solution for low risk users.

### Key features

- **Ideal for low care applications** - suited where support and comfort is required and individuals are at low risk of pressure injury
- **High stability** - quality foam provides an effective support surface for patient positioning and transferring
- **Advantiflex cover** - high-performance, single side PU waterproof cover. Cover is a 4-way stretch material for shear protection, immersion and envelopment
- **Construction** - all seams high frequency sealed and welded, waterfall fluid penetration flap
- **Integrated zipper** - allows easy removal of cover, robust and durable

## ADDITIONAL INFORMATION

<b>Sizes</b>	Single, King Single
<b>Safe Working Load (SWL)</b>	150 kg
<b>Length</b>	198 cm
<b>Width</b>	88 cm
<b>Depth</b>	15 cm
<b>Brand</b>	Forte Healthcare