



SKU: S1040

DRIVE MEDICAL STELLA ELECTRIC RECLINER LIFT CHAIR - LEATHER - TWIN MOTOR

\$2,499.99

The Days Healthcare Stella Lift Chair is an electrically operated, fully upholstered lounge chair with dual motors. Comes equipped with a **2 year comprehensive warranty**. The Stella Twin Motor Riser Recliner Lift Chair is a reliable recliner with great additional features including a seat with foam padding and a two section padded backrest for lumbar and head support.

- Two-section padded backrest offering lumbar support and a comfortable headrest.
- Swivel castors for easy movement
- Dual-Motor chair features two motors.
- It will rise and tilt users to their feet.
- Quiet motors provide a smooth and discreet operation.
- Ergonomically contoured seat and back.
- Fully padded chaise for added comfort.
- Two convenient side pockets for easy storage.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	s2051	\$2,499.99	Navy (Light Charcoal as pictured) Leather	
	S1040	\$2,499.99	Chocolate Leather	

PRODUCT DESCRIPTION

- Quality assured, meeting all relevant fire safety standards.
- Built-in fail-safe battery back-up system.*
- Comes with a simple-to-operate handset.
- Swivel castors enable easy manoeuvrability of the chair.
- Stitched handset holder to prevent loss / damage to handset.
- Ideal for individuals who have difficulty standing up / sitting down.
- 2-year warranty on wooden frame and electrical components.
- Maximum weight capacity 130kg (20st / 285lb).

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Brand	Drive Medical
Seat Height	470 mm
Seat Depth	500 mm
Seat Width	530 mm
Overall Width	950 mm
Overall Height	1200 mm
Max SWL	130 kg
Material	Leather
Wattage	60W