



**SKU:** S1040

## **DRIVE MEDICAL STELLA ELECTRIC RECLINER LIFT CHAIR - LEATHER - TWIN MOTOR**

**\$2,499.99**



The Days Healthcare Stella Lift Chair is an electrically operated, fully upholstered lounge chair with dual motors. Comes equipped with a **2 year comprehensive warranty**. The Stella Twin Motor Riser Recliner Lift Chair is a reliable recliner with great additional features including a seat with foam padding and a two section padded backrest for lumbar and head support.

- Two-section padded backrest offering lumbar support and a comfortable headrest.
- Swivel castors for easy movement
- Dual-Motor chair features two motors.
- It will rise and tilt users to their feet.
- Quiet motors provide a smooth and discreet operation.
- Ergonomically contoured seat and back.
- Fully padded chaise for added comfort.
- Two convenient side pockets for easy storage.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Colours
	s2051	\$2,499.99	Navy (Light Charcoal as pictured) Leather	
	S1040	\$2,499.99	Chocolate Leather	

## PRODUCT DESCRIPTION

- Quality assured, meeting all relevant fire safety standards.
- Built-in fail-safe battery back-up system.\*
- Comes with a simple-to-operate handset.
- Swivel castors enable easy manoeuvrability of the chair.
- Stitched handset holder to prevent loss / damage to handset.
- Ideal for individuals who have difficulty standing up / sitting down.
- 2-year warranty on wooden frame and electrical components.
- Maximum weight capacity 130kg (20st / 285lb).

### **Single vs twin vs quad motor lift chairs - what's the difference?**

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## ADDITIONAL INFORMATION

<b>Brand</b>	Drive Medical
<b>Seat Height</b>	470 mm
<b>Seat Depth</b>	500 mm
<b>Seat Width</b>	530 mm
<b>Overall Width</b>	950 mm
<b>Overall Height</b>	1200 mm
<b>Max SWL</b>	130 kg
<b>Material</b>	Leather
<b>Wattage</b>	60W