



**SKU:** S1111

## SUCTION GRAB RAIL

**\$69.99**

This Suction Grab Rail is great for providing extra support in the toilet, shower and bathroom. It is easy to install and can be easily removed.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

Simply press down on one lever and then the other to secure in place. No drilling necessary. Ideal for using at home or to take traveling.

## ADDITIONAL INFORMATION

**daily-living**

Bathroom Other

**Handgrip Length**

500 mm or 600 mm

**Brand**

ILS ADL