



SKU: S1111

HEMOCRAFT SUCTION GRAB RAIL

\$69.99

The Homecraft Suction Grab Rail provides tool-free, flexible support for safer movement in bathrooms and similar areas. Its strong suction cups and easy-lock system make it a reliable, non-permanent aid on smooth, nonporous surfaces like tiles or glass.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Homecraft Suction Grab Rail enhances safety in the bathroom by providing a secure, temporary handhold. Designed for use on smooth, nonporous surfaces such as tiles, glass, porcelain or fibreglass, it attaches easily using large suction pads – no screws or drilling required.

This rail features quick-release levers and safety lock indicators that change colour if the suction weakens, offering peace of mind during use.

The textured surface offers a stable grip even with wet hands, while the lightweight, portable design makes it easy to reposition or take on the go.

It's an ideal solution for short-term use, travel or rental homes where permanent fixtures are not an option.

Who it's for

- For clients seeking temporary support during bathroom activities.
- For therapists recommending non-permanent grab rails to enhance client safety and independence in various settings.

Key features

- **Tool-free suction installation** – attaches securely to smooth tiles or glass
- **Safety lock indicators** – visually signal when suction strength is reduced
- **Textured grip surface** – enhances grip in wet bathroom environments
- **Quick-release levers** – allow easy repositioning or removal
- **Lightweight and portable** – perfect for travel, temporary use or rental homes
- **No drilling required** – preserves bathroom surfaces

ADDITIONAL INFORMATION

Maximum Weight Capacity	63 kg (Straight/Angled Pull)
Length	50 cm