

**SKU:** S11755

## **SUNRISE MEDICAL STERLING S425 ELECTRIC MOBILITY SCOOTER**

**\$3,999.99**

For those who fancy a Fast, strong and comfortable long ride, then the Sunrise Medical Sterling S425 mobility scooter is the one for you,

The larger battery size of 50AH will help to power a journey up to 42 kms, packed with brilliant features and superb driving dynamics such as all-round suspension and large 30 cm pneumatic wheels. Combining bright, front/rear advanced LED lights with a maximum speed of 12km/h, the S425 mobility scooter is also approved for road use



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

### **Sunrise Medical Sterling S425 Electric Mobility Scooter**

The Sterling S425 mobility scooter has been designed with advanced seating solutions and intuitive controls so you experience comfort on every trip. It is perfect for indoor or outdoor use - day or night! Turn heads with the S425's unique carbon-metallic colour just as easily as the S425 mobile scooter does.

#### **Key features include:**

- All-round suspension – The larger 30 cm pneumatic wheels and all-around adjustable suspension system ensures you experience a comfortable ride when navigating over bumpy terrain or tackling those tricky outdoor obstacles. To achieve an even firmer ride, there is an optional heavy-duty rear-suspension kit for the S425 mobility scooter.
- Advanced LED lighting technology - the lights on the S425 mobility scooter are 400 times more energy efficient than standard bulbs so you use less battery power on every ride.
- Easy-to-adjust luxury seating - Adjust the seat height, depth and recline to suit your own unique driving position.
- Flip-up comfort armrests which are width, angle and depth adjustable.
- Seat rotation is standard to enable an effortless transfer in and out of the S425 mobility scooter.
- Easy-to-reach ergonomic delta handle-bar.
- Adjustable tiller – the S425 mobility scooter's clever design even ensures that the headlight is automatically tilted as the tiller is adjusted.
- Unique wig-wag control system using your fingers or thumb.
- Waterproof soft touch controls such as easy-to-use speed adjustment, LED battery gauge and a powerful horn.
- Removable basket with carry handle.
- 2 carbon effort mirrors.

Additional accessories which can be purchased separately include:

- Canopy
- Rear bag and cane holder
- Rear trailer

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	150 kg (330.69 lbs)
<b>Product Battery Size / Type</b>	2 x 50 AH Gel
<b>Product Range</b>	42 km
<b>Product Speed (Max)</b>	10 km/hr
<b>Product Colour</b>	Carbon metallic
<b>Product Warranty</b>	Mechanical 24 months Electrical 24 months Batteries 6 months
<b>Product Size - in carton(s) (WxLxH) in cm</b>	72 x 144 x 87
<b>Product Size - in carton(s) (gross weight) in kg</b>	150 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	66 x 135 x 129
<b>Product Size - overall assembled (gross weight) in kg</b>	117 kg
<b>Product Size - seat base width in cm</b>	49 cm
<b>Product Size - seat base depth in cm</b>	49 cm

<b>Product Size - seat back width in cm</b>	47.5 cm
<b>Product Size - seat back height in cm</b>	52.5 cm
<b>Product Size - seat to floor in cm</b>	45 cm – 50 cm
<b>Product Size - front/ rear tyre type/size</b>	11.8" / Pneumatic
<b>Product Size - front/rear wheel size in inches</b>	5"
<b>Product Drive - turning circle in cm</b>	200 cm
<b>Product Drive - max safe ground clearance in cm</b>	6 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	Yes (Front and Rear)
<b>Product Charger Type</b>	8 AH
<b>Product Motor Size</b>	4 Pole / 1200 watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Analog