

SKU: S11755

SUNRISE MEDICAL STERLING S425 MOBILITY SCOOTER

\$4,299.99

The Sunrise Medical Sterling S425 Mobility Scooter is a versatile solution for indoor and outdoor use. With all-round suspension, 12" wheels and energy-efficient lighting, it balances comfort, control and convenience.



Scan QR code to
view product online

GALLERY IMAGES





PRODUCT DESCRIPTION

The Sterling S425 is designed to offer reliable performance across everyday settings. It features 12" pneumatic wheels and all-round suspension to manage rough terrain.

Its luxury vinyl seat includes depth and recline adjustments, as well as a swivel base and adjustable armrests and headrest. The unique wig-wag throttle allows control using fingers or thumbs, ideal for those with limited hand function.

It features a bright, efficient LED lighting system, ergonomic tiller and stylish metallic silver finish.

The Sterling S425 includes a seat belt, front basket and two rear-view mirrors as standard.

Who it's for

- For clients needing a comfortable scooter for indoor and outdoor use.
- For therapists recommending a versatile scooter with accessible controls and ergonomic support.

Key features

- **All-round suspension with 12" pneumatic wheels** - supports comfort over uneven surfaces
 - **Advanced LED lighting system** - increases safety while conserving battery
 - **Luxury swivel seat with multiple adjustments** - supports comfort and transfers
- **Wig-wag control system and delta tiller** - easy handling for those with reduced hand function
- **Stylish metallic silver finish with seat belt and mirrors** - safety and style combined

ADDITIONAL INFORMATION

Safe Working Load (SWL)	150 kg
Battery Type	Gel
Battery Size	2 x 50 AH
Range	42 km
Max Speed	10 km/h
Size (W x L x H assembled)	66 cm x 135 cm x 129 cm
Weight	117 kg
Suspension	Front & Rear
Display	Analog
Max Safe Incline Slope	8°
Turning Circle	200 cm
Ground Clearance	6 cm
Brand	Sunrise Medical

[Sunrise Medical Sterling S425 Brochure](#)

[Sunrise Medical Sterling S425 User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

WHAT IS THE TOTAL WIDTH OF THE SUNRISE MEDICAL S425 MOBILITY SCOOTER?

The overall width for the Sunrise Medical S425 Mobility Scooter is 660mm (66cm)

DOES THIS GET DELIVERED ASSEMBLED AND ARE THE WHEELS LARGER FOR GOING SAFELY UP AND DOWN CURBS

S425 Mobility scooter needs minimal assembly. All you have to do is add Seat to the scooter. S425 Mobility scooter has 12" pneumatic tyres.