

**SKU:** S11755

## **SUNRISE MEDICAL STERLING S425 ELECTRIC MOBILITY SCOOTER**

**\$4,299.99**

For those who fancy a Fast, strong and comfortable long ride, then the Sunrise Medical Sterling S425 mobility scooter is the one for you,

The larger battery size of 50AH will help to power a journey up to 42 kms, packed with brilliant features and superb driving dynamics such as all-round suspension and large 30 cm pneumatic wheels. Combining bright, front/rear advanced LED lights with a maximum speed of 12km/h, the S425 mobility scooter is also approved for road use



Scan QR code to  
view product online



## PRODUCT DESCRIPTION

### Sunrise Medical Sterling S425 Electric Mobility Scooter

The Sterling S425 mobility scooter has been designed with advanced seating solutions and intuitive controls so you experience comfort on every trip. It is perfect for indoor or outdoor use - day or night! Turn heads with the S425's unique carbon-metallic colour just as easily as the S425 mobile scooter does.

#### Key features include:

- All-round suspension – The larger 30 cm pneumatic wheels and all-around adjustable suspension system ensures you experience a comfortable ride when navigating over bumpy terrain or tackling those tricky outdoor obstacles. To achieve an even firmer ride, there is an optional heavy-duty rear-suspension kit for the S425 mobility scooter.
- Advanced LED lighting technology - the lights on the S425 mobility scooter are 400 times more energy efficient than standard bulbs so you use less battery power on every ride.
- Easy-to-adjust luxury seating - Adjust the seat height, depth and recline to suit your own unique driving position.
- Flip-up comfort armrests which are width, angle and depth adjustable.
- Seat rotation is standard to enable an effortless transfer in and out of the S425 mobility scooter.
- Easy-to-reach ergonomic delta handle-bar.
- Adjustable tiller – the S425 mobility scooter's clever design even ensures that the headlight is automatically tilted as the tiller is adjusted.
- Unique wig-wag control system using your fingers or thumb.
- Waterproof soft touch controls such as easy-to-use speed adjustment, LED battery gauge and a powerful horn.
- Removable basket with carry handle.
- 2 carbon effort mirrors.

Additional accessories which can be purchased separately include:

- Canopy
- Rear bag and cane holder
- Rear trailer

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	150 kg (330.69 lbs)
<b>Product Battery Size / Type</b>	2 x 50 AH Gel
<b>Product Range</b>	42 km
<b>Product Speed (Max)</b>	10 km/hr
<b>Product Colour</b>	Carbon metallic
<b>Product Warranty</b>	Mechanical 24 months Electrical 24 months Batteries 6 months
<b>Product Size - in carton(s) (WxLxH) in cm</b>	72 x 144 x 87
<b>Product Size - in carton(s) (gross weight) in kg</b>	150 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	66 x 135 x 129
<b>Product Size - overall assembled (gross weight) in kg</b>	117 kg
<b>Product Size - seat base width in cm</b>	49 cm
<b>Product Size - seat base depth in cm</b>	49 cm

<b>Product Size - seat back width in cm</b>	47.5 cm
<b>Product Size - seat back height in cm</b>	52.5 cm
<b>Product Size - seat to floor in cm</b>	45 cm – 50 cm
<b>Product Size - front/ rear tyre type/size</b>	11.8" / Pneumatic
<b>Product Size - front/rear wheel size in inches</b>	5"
<b>Product Drive - turning circle in cm</b>	200 cm
<b>Product Drive - max safe ground clearance in cm</b>	6 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	Yes (Front and Rear)
<b>Product Charger Type</b>	8 AH
<b>Product Motor Size</b>	4 Pole / 1200 watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Analog