

SKU:

SUNRISE MEDICAL STERLING S700 MOBILITY SCOOTER

\$4,999.99

The S700 mobility scooter is the leader of the pack with superior outdoor performance features. Tackle greater gradients of 10 degrees, higher kerbs of up to 100 mm and greater obstacles with larger 33 cm tyres – it's truly a mobility scooter designed for those who just can't sit still. There's no need to worry about distance either as the large 75Ah battery size will power an extended range of up to 55 km.

(2x75AH MK Gel Batteries)



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Product Colour
	S11756	\$4,999.99		Silver

PRODUCT DESCRIPTION

Pure outdoor performance.

Experience a smooth and rewarding drive each and every time with the S700's all-round advanced suspension. Every wheel has its own independent suspension, so whether you're cruising around town at 6 kph or travelling over more varied terrain at 10 kph, you're assured of the ultimate comfort over every bump. Plus, with bigger 33 cm pneumatic wheels, you'll effortlessly clear any obstacles with more confidence!

Powerful LED's to light your way.

With advanced LED lighting technology, the lights on the S700 mobility scooter are 400 times more energy efficient than standard bulbs, so not only are you using less battery power on every trip, you also won't have to worry about changing bulbs either.

Comfort on every trip.

With 75 Ah high performance batteries powering a potential range of 55 km, finding your perfect seating position for those long journey's is effortless with the easy-to-adjust luxury seating. Adjust the seat height, depth and recline to suit your own unique driving position. The flip-up comfort armrests are also width, angle and depth adjustable - and combined with the seat rotation ensures that transferring in and out of your mobility scooter is effortless.

Intuitive controls

Seeing as you'll be spending a lot of time behind the controls of your S700, we've ensured the delta handle-bar is ergonomic and easily within reach thanks to the infinitely adjustable tiller (the clever design even ensures that the headlight is automatically tilted as the tiller is adjusted). The unique wig-wag control system also allows you to control your mobility scooter by using your fingers or thumb

ADDITIONAL INFORMATION

Product Safe Working Load (SWL)	160 kg (352.7 lbs)
Product Battery Size / Type	2 x 75 AH Gel
Product Range	55 km
Product Speed (Max)	10 km/hr
Product Colour	Silver
Product Warranty	Mechanical 24 months Electrical 24 months Batteries 6 months
Product Size - in carton(s) (WxLxH) in cm	74 x 154 x 90
Product Size - in carton(s) (gross weight) in kg	160 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	67 x 150 x 129
Product Size - overall assembled (gross weight) in kg	146 kg
Product Size - seat base width in cm	49 cm
Product Size - seat back width in cm	47.5 cm

Product Size - seat back height in cm	52.5 cm
Product Size - seat to floor in cm	45 cm – 50 cm
Product Size - front/ rear tyre type/size	13" / Pneumatic
Product Size - front/rear wheel size in inches	5"
Product Drive - turning circle in cm	220 cm
Product Drive - max safe ground clearance in cm	10 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	Yes (Front and Rear)
Product Charger Type	8 AH
Product Motor Size	4 Pole / 2700 watt
Product Lighting	LED
Product Display	Analog