



SKU:

SUNRISE MEDICAL XPLORE 2

Rear-Wheel Powered Wheelchair All Terrain Hybrid Design

The Quickie® Xplore 2™ is the most hybrid product available. This wheelchair offers the best of both worlds; the traction of a rear wheel drive powered wheelchair and the short turning radius of a mid-wheel drive wheelchair.

The fast speed of the Xplore 2™ makes it the ideal choice for someone spending a lot of time outdoors and who wants to go the distance.

Safe and comfortable, this motorised wheelchair feature a 6-wheel independent suspension which is available in both the Standard and Heavy Duty models. The electric wheelchair come standard with transport brackets and has a wide range of seating options to choose from - tilt, manual & power recline, and elevating seat.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Xplore is an excellent choice for an active user spending time outdoors. The directional stability at high speeds in all types of environments, helps ensure that the end user's optimal Center of Gravity will be maintained without any adjustments required.

ADDITIONAL INFORMATION

Weight Capacity	136 KG Standard, 181 KG Heavy Duty option
Front Caster	8" Caster, 9" Caster pneumatic with or without insert
Drive Wheel	14" Cast wheel pneumatic with or without insert
Transit	c4-Point wheelchair transport included
Electronics	R-Net EL 90 Amp fully expandable controller R-net 120 Amp fully expandable controller OPTION
Seat to floor height	16.5"
Max Speed	14 KPH
Overall Base Width	65 CM
Overall Base Length	110 CM
Seating	Fully adjustable seat: Eclipse
Warranty	5 years for frame, please see owner's manual for all warranty information
Brand	Sunrise Medical