



**SKU:**

## **INVACARE HYGIENE LOW FABRIC SLING**




**\$399.99**

The Invacare Hygiene Low Fabric Sling offers safe, simplified toileting transfers for individuals with adequate trunk, hip and head control. Its wide opening and low-back design support comfort, access and ease of use.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S306434	\$399.99		X-Large
	S300785	\$399.99		Medium
	S300784	\$399.99		Large

## PRODUCT DESCRIPTION

The Invacare Hygiene Low Fabric Sling is a specialised toileting sling designed for seat-to-seat transfers. Made from durable polyester, it provides a firm yet comfortable surface that slides easily into place and dries quickly after washing – ideal for regular hygiene routines.

The sling's wide toileting aperture allows clothing to be removed easily while maintaining support around the lower back and hips. It is intended for those who have sufficient control of their trunk, hips and head.

Features include two central alignment markers for correct spinal positioning, colour-coded loops for safe and accurate attachment, and a supportive body belt for additional security.

A backrest positioning pocket and carefully placed handles assist carers during fitting and repositioning.

### Who it's for

- For clients with good trunk, hip and head control requiring secure, seated toileting transfers.
- For therapists seeking a toileting sling with low-back support and quick-application design.

### Key features

- **Low-back design** – provides support for those with independent head and trunk control
- **Wide toileting aperture** – allows clothing access for hygiene routines
- **Durable polyester fabric** – slides easily, dries quickly and withstands frequent washing
- **Positioning handles and backrest pocket** – assist carers during application
- **Centre line markers** – help align spine for safe lifting
- **Coloured loop options** – support accurate and secure attachment
- **Adjustable belt** – provides added support and positioning control
- **Machine washable up to 90°C** – supports hygiene protocols
- **Available in 3 sizes** – with a safe working load up to 250 kg

## ADDITIONAL INFORMATION

<b>Brand</b>	Invacare
<b>Size</b>	Medium, Large, X-Large
<b>Material</b>	Fabric
<b>Sling Type</b>	Hygiene