


SKU: S1798

DAYS HEALTHCARE SWIFT 18" WHEELCHAIR

\$549.99

The Swift Transit Wheelchair has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional chair for both indoor and outdoor use. Its compact dimensions allow good maneuverability indoors and the ultra light construction and ease to fold for storage or transport make it a very versatile chair.



Scan QR code to
view product online



PRODUCT DESCRIPTION

With an ultra lightweight aluminium frame, which easily folds and dismantles without any tools, the Swift wheelchair is perfect for occasional use. Designed for days in or out, it's ideal for taking on trips, leaving at relatives' homes, or even storing in the cupboard so it's there in case you need it.

Features

- Ultra light frame 7.8kg
- Folding backrest (for storage or transport)
- Quick release wheels (for transport)
- Fixed arms
- Swing away and removable leg rests
- Solid tyres (no maintenance required)
- Extended brake handle (ease of brake application)

ADDITIONAL INFORMATION

Brand	Days Healthcare
Max SWL	115 kg
Frame material	Aluminium
Seat height - front	500 mm
Seat height - rear	475 mm
Seat width	460mm
Seat Depth	420 mm
Seat angle	3.5°
Backrest height	410 mm
Backrest angle	8°
Legrest angle	70°
Armrest height	230 mm
Chair weight	7.8kg (frame only) Weight: 12.9kg (complete)