

SKU:

TALLEY PULSAIR CHOICE MATTRESS

~~\$5,999.99~~ \$5,999.99



The Pulsair Choice range is designed to bring maximum flexibility of pressure area care to patients in the home care environment. Mix and match the PULSAIR CHOICE power unit with the mattress replacement for simple, cost-effective stepping up or down of pressure care for patients at all risks of pressure damage.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S6927	\$4,499.99		Single Standard
	S9714	\$4,499.99		Single Long
	S311258	\$5,999.99		King Single Standard
	S311570	\$5,999.99		King Single Long

PRODUCT DESCRIPTION

Pulsair Choice products provide a clinical preventative measure to provide pressure relief and manage patients at risk of developing or with existing pressure Injuries. Please see below some of the product's key features:

- **ODSTM cell design for optimum patient support and security:** Air cells feature ortho differential support (ODSTM) which provides optimum internal cell pressures for varying body weights and shapes, automatically compensating for the patient's BMI. Narrower in the centre than at the outer edge, air cells feature in-built support pressure differential, supporting the smaller, lighter patient on the softer central area of the cell, and the larger, heavier patient on the firmer outer edges, as well as providing a feeling of security for the patient
- **AUTO support - Fully automatic pressure adjustment:** The AUTO support feature, automatically adjusts pressure within each of the cells throughout the cycle so that support, posture, and therapy are constantly maintained at optimum levels, in response to patient weight, movement and position. Equalisation of cell pressure automatically takes place at each stage of the cycle, again to ensure precise pressure and therapy is applied.
- **Comfort Control:** Allows the carer to adjust the comfort setting in therapy mode to help patient compliance.
- **Easy to Operate:** Simple, concise operation buttons with auto lock feature to prevent inadvertent button operation. Power unit features audio alerts together with LCD screen showing current operation status / fault conditions.
- Ideal for home care environments
- Removable/full launderable cover
- CPR facility
- Patient transport facility
- Vapour-permeable, waterproof, two way stretch cover

ADDITIONAL INFORMATION

Brand	Direct Healthcare Group
Length	1950 mm
Width	890 mm
Depth	180 mm
Inner Construction Materials	ODS deep cell PU air cells set
Cover Material	PU coated stretch nylon
Number of Cells	30
Max. User Weight	200 kg
Min. User Weight	None
SWL	200kg