

**SKU:**

## **TILITE AERO T**



The Aero T is a highly adjustable and rigid aluminium wheelchair that boasts increased push efficiency and a stylish, horizontal bottom tube for easy car transfers.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

If you're in the market for a highly adjustable and rigid aluminium wheelchair, the Aero T from TiLite is an excellent option. With its dual-tube geometry, the Aero T provides increased push efficiency, which means you'll need fewer push strokes throughout the day. Additionally, the Aero T features a stylish, horizontal bottom tube that allows for easy car transfers. But perhaps the best feature of the Aero T is its full adjustability in both the Centre of Gravity and Rear/Front Seat to Floor Height, allowing you to customise the chair to your specific needs and preferences.

## ADDITIONAL INFORMATION

<b>Brand</b>	TiLite
<b>Frame Type</b>	Rigid
<b>Material</b>	Aluminum
<b>Weight Capacity</b>	120 kg (Standard), 159 kg (Heavy Duty)
<b>Transit Approved</b>	No
<b>Transport Weight</b>	5.4 kg
<b>Pediatric Seating</b>	No
<b>Seat Width Range</b>	12" - 20"
<b>Seat Depth Range</b>	12" - 20"